

Overview

► **Purpose** – To determine if MFLC (Microflow Liquid Chromatography) coupled with an MS system can generate bioanalytical validation data that exceeds FDA criteria and is comparable to conventional HPLC-MS/MS systems.

► **Methods** – Validation of Methotrexate in human plasma utilizing MFLC-MS/MS (Eksigent ekspert™ microLC 200 System coupled with an ABSciex 5500® QTRAP) and HPLC-MS/MS (Shimadzu LC-10AD coupled with an ABSciex 5500® QTRAP).

► **Results** – The MFLC-MS/MS system provides bioanalytical validation data that passes FDA criteria and is comparable to the HPLC-MS/MS validation data.

Introduction

Micro-flow liquid chromatography coupled with a mass spectrometer (MFLC-MS/MS) has many advantages (sensitivity increase, less solvent consumption and reduced source contamination). These advantages have led to the manufacture of specialized LC pumps that accurately deliver ≤ 100 $\mu\text{L}/\text{min}$ of solvent and LC columns with an inside diameter of ≤ 1 mm. In order to determine if the MFLC-MS/MS system can provide bioanalytical data that exceeds the criteria required by the FDA, a validation (determination of Methotrexate in human plasma) was performed. A validation was also performed using the conventional HPLC-MS/MS for comparison purposes.

Here we report data generated for the determination of Methotrexate from human plasma by MFLC-MS/MS and HPLC-MS/MS.

Methods

Extraction

- 50 μL of human plasma sample
- 25 μL of the internal standard Methotrexate- D_3
- Acetonitrile added for protein precipitation
- Supernatant transferred and water 1% formic acid added

MFLC

- Eksigent ekspert™ microLC 200 System gradient MFLC using acetonitrile and water with 1% formic acid
- Flow rate: 35 $\mu\text{L}/\text{minute}$
- Column: ProntoSIL 120-3-C18, 0.5 x 50 mm (MAC-MOD)
- Column temperature: 50°C

Conventional HPLC

- Shimadzu LC-10AD using acetonitrile and water with 1% formic acid
- Flow rate: 700 $\mu\text{L}/\text{minute}$
- Column: ProntoSIL 120-3-C18, 2.0 x 50 mm (MAC-MOD)
- Column temperature: 50°C

MS

- ABSciex 5500® QTRAP operating in MRM mode
- ESI
- Positive ion mode
- MRM transitions:
 - Methotrexate: 455.4 \rightarrow 308.1
 - Methotrexate- D_3 : 458.4 \rightarrow 311.5





Table 1: Validation Data Comparison of HPLC and MFLC Systems

	Shimadzu HPLC 700 $\mu\text{L}/\text{min}$	Eksigent MFLC 35 $\mu\text{L}/\text{min}$
*Intraday A/P (%)	107 \pm 11.8	93.9 \pm 12.9
*Interday A/P (%)	112 \pm 13.8	106 \pm 13.4
Matrix Factor	1.02	1.07

* All data includes the LLOQ and ULOQ values




Figure 1. MFLC-MS/MS and HPLC-MS/MS chromatograms from the analysis of a 2.0 μL injection of Methotrexate extracted from human plasma.

Conclusions

- The MFLC-MS/MS system can be used to perform bioanalytical validations.
- The validation data produced by the MFLC-MS/MS system is equivalent to the data generated from the conventional HPLC-MS/MS system.
- Only ~ 125 mL of solvent was needed to perform the MFLC-MS/MS validation compared to $\sim 2,500$ mL of solvent required for the HPLC-MS/MS validation.
- MFLC-MS/MS system is less susceptible to contamination than the HPLC-MS/MS system.




Figure 2. QTRAP® 5500 interface plate after ~ 400 extracted human plasma injections from the Eksigent system – no switching valve.




Figure 3. QTRAP® 5500 interface plate after ~ 150 extracted human plasma injections from the Shimadzu conventional HPLC system – no switching valve.



Figure 4. AB SCIEX QTRAP® 5500 Coupled with an Eksigent ekspert™ microLC 200 MFLC System.