

OVERVIEW

Purpose

Non-invasive matrices such as tears, sweat, saliva, and milk offer an enticing alternative to traditional sampling from blood, serum, and plasma. Collection of these matrices is simpler and more affordable than venipuncture. It is also much less unpleasant for patients and study subjects allowing for more frequent sampling, higher compliance, and removal of fear of needles as a barrier to clinical trial recruitment. However, bioanalysis of these matrices is not without challenges. Sample volumes can be very low, the more complex matrices such as milk may require extensive sample preparation, and the use of the matrix must be biologically relevant. Using the example of a method we developed at Alturas Analytics Inc. to measure Tobramycin from human tears, we will discuss practical applications of microsampling of non-invasive matrices for bioanalysis, and improvements offered by emerging technologies.

Objectives

- Understand the benefits and challenges of using non-invasive matrices in bioanalysis
- Explore currently available microsampling technologies and advances in analytical techniques
- Discuss strategies for collecting, extracting, and analyzing non-invasive matrices

Microsampling of Non-Invasive Matrices: Practical Examples Using Tears and a Perspective of Past and Emerging Technologies

Presented by Jason Watts, Ph.D.

Alturas Analytics, Inc., Moscow, ID, United States

INTRODUCTION

Microsampling

- » Typically uses $\leq 50 \mu\text{L}$
- » Simple collection and storage
- » Possibility for at-home sampling
- » Site-Centric \rightarrow Patient-Centric


Non-invasive Matrices

- » Tears, sweat, saliva, milk
- » High correlation with non-protein bound plasma concentrations for many drugs


Practical Considerations

- » Can sufficient sensitivity be achieved?
- » Analyte stability
- » Does the matrix concentration correlate with plasma concentrations?
- » What is gained over traditional sampling?

Adapting blood microsampling techniques for non-invasive matrices



- » Improved precision and accuracy using color-indicating dyes



- » Improving extraction workflow with common Bioanalytical tools

MATRICES

Tears

- » Constantly produced at $1.2 \mu\text{L}/\text{min}$
- » Unstimulated volume $\sim 7 - 10 \mu\text{L}$
- » 0.6 - 0.8% protein (0.4% albumin)
- » pH 6.5 - 7.6
- » Differences between stimulated and unstimulated collection





Capillary tube collection
From Posa et al., 2012



Schirmer strip collection
From Posa et al., 2012

Tobramycin from Human Tears



QC Level ($\mu\text{g}/\text{mL}$)	Assay Accuracy and Precision ($\% \pm \%CV$)	Matrix Factor
12	90.0 ± 4.5	NA
3.0	96.2 ± 0.07	NA
1.5	101 ± 5.7	0.97

CONTINUED FROM FRONT

Sweat

- » Collected on sweat wipes or patches
- » Volume variability: normalized to sodium and potassium levels
- » pH 4 - 6.8 when resting
 - High sweat/plasma ratio for basic drugs
- » Commonly used for monitoring drugs of abuse

Saliva

- » Production and composition
 - 0.5 mL/min
 - pH ~ 6 - 7
 - Excretion of drugs dependent on permeability and protein binding
- » Production and composition
 - Unstimulated
 - Stimulated

Using blood microsampling device for saliva sampling

- » Mitra® (Neoteryx) microsampling device
 - Precise sample collection
 - Ideal for pediatric or at home sampling
 - 10, 20 or 30 µL sample volume





Fig. 2. Correlation between concentrations of dexamethasone in plasma and saliva taken simultaneously from healthy subjects and patients pretreated with dexamethasone.


Correlation between plasma and saliva dexamethasone concentrations
Thijssen et al., 1996



SIMPLE EXTRACTION OF DEXAMETHASONE FROM SALIVA USING MITRA® MICROSAMPLING DEVICE


Method

- » Collect sample on 20 µL Mitra® tip
- » Dry for at least 2 hours at ambient temperature
- » Place dried tip in 96 well plate with curve and QCs
- » Add 25 µL 100 ng/mL IS (Dex-D4 in ACN:H2O 1:1)
- » Add 300 µL ACN. Vortex and incubate 30 min
- » Transfer extract to fresh plate. Evaporate to dryness
- » Reconstitute in 100 µL of ACN:H2O with 0.1% formic acid



0.4 ng/mL Dexamethasone extracted from Saliva

Dexamethasone extracted from saliva



Dexamethasone 392.5 g/mol
Structure from <https://pubchem.ncbi.nlm.nih.gov/compound/Dexamethasone>

Concentration (ng/mL)	Assay Accuracy and Precision (% ± % CV)	% Recovery
0.8	98.3 ± 0.94%	94
50	103.5 ± 0.68%	117
400	100.1 ± 7.0%	121

REFERENCES

- Francesca G. Bellagambi, Tommaso Lomonaco, Pietro Salvo, Federico Vivaldi, Marie Hangouet, Silvia Ghimenti, Denise Biagini, Fabio Di Francesco, Roger Fuoco, And Abdelhamid Errachid., Saliva Sampling: Methods and devices. An overview. Trends in Analytical Chemistry 124 115781. 2020
- Andreas Posa, Lars Brauer, Martin Schicht, Fabian Garreis, Stephanie Beileke, and Friedrich Paulsen., Schirmer strip vs. capillary tube method: Non-invasive methods of obtaining proteins from tear fluid. Annals of Anatomy 195 137-142. 2013
- Jos H.H. Thijssen, Christine C.Gispens-de Wied, Gemma M. Van Heeswijk, and Winnie Veeman., Determination of dexamethasone in saliva. Clinical Chemistry 42:8 1238-1242. 1996



Alturas Analytics is a GLP compliant bioanalytical CRO specializing in MS/MS solutions supporting early discovery through late phase clinical trials.

In addition to providing PK support services to pharmaceutical companies worldwide, Alturas maintains an intensive research effort of applying new technologies leading to scientific advancement.